

PRODIGY MEDICAL

Phone: (907) 226-3400 - Fax: (907) 226-3300 - Email: admin@prodigyalaska.com
 347 West Danview Avenue – PO Box 2624 – Homer AK 99603

SIMPLYGO BATTERY LIFE CHART

		PULSE DOSE SETTING					
		1	2	3	4	5	6
NUMBER OF BATTERIES	1	3.7 hours	3.5 hours	3.1 hours	2.7 hours	2.3 hours	1.9 hours
	2	7.4 hours	7 hours	6.2 hours	5.4 hours	4.6 hours	3.8 hours
	3	11.1 hours	10.5 hours	9.3 hours	8.1 hours	6.9 hours	5.7 hours
	4	14.8 hours	14 hours	12.4 hours	10.8 hours	9.2 hours	7.6 hours

		CONTINUOUS FLOW SETTING		
		0.5	1	2
NUMBER OF BATTERIES	1	3.1 hours	2.3 hours	0.7 hours
	2	6.2 hours	4.6 hours	1.4 hours
	3	9.3 hours	6.9 hours	2.1 hours
	4	12.4 hours	9.2 hours	2.8 hours